

Weasis for Windows Guide

VERSION 1.2

ZED TECHNOLOGIES PTY. LTD. 2017

TABLE OF CONTENTS

Disclaimers 3

History 4

Extracting study & opening Weasis in windows 5

DISCLAIMERS

EXAMPLE IMAGES USED

Images contained within this documentation are copyright © Zed Technologies Pty. Ltd.

HISTORY


The history of this document is outlined below. Zed Technologies strives to ensure the documentation is kept up to date.

Version	Author	Date
V1.2	Oscar Putra	29/05/2017

EXTRACTING STUDY & OPENING WEASIS IN WINDOWS

Once export study download has been completed, the user must locate the file, press right-click on the file to unzip/extract it as it's currently in zip format and cannot be used until it has been unzipped or extracted.

After file has been extracted, open the folder and double click **viewer-win32.exe** to open DICOM medical viewer and view the exported study.



In Weasis for Windows, the user can double click on the left pane to select which series to view and use the mouse scroll wheel to scroll through images.

